



100bp DNA Ladder plus (ready to use)-100-3000bp

Cat. No. SL7041

Size: 50µg / 500µl

Description:

A unique combination of PCR products and a number of proprietary plasmids digested with appropriate restriction enzymes to yield 12 fragments, suitable for use as molecular weight standards for agarose gel electrophoresis. The DNA includes fragments ranging from 100-3,000 base pairs. The 500 and 1,500 base pair bands have increased intensity to serve as reference points. The approximate mass of DNA in each band is provided (0.54 µg a load) for approximating the mass of DNA in comparably intense samples of similar size.

Source:

PCR products and double-stranded DNA digested with appropriate restriction enzymes, are phenol extracted and equilibrated to 10 mM Tris-HCl (pH 8.0) and 10 mM EDTA.

Range: 100-3,000 bp

Number of bands: 12

Concentration: 108 µg/ml

Package: 50µg / 500µl

Recommended Load: 5 µl / well


Containing orange G & xylene cyanol FF as tracking dyes.

Storage:

Store at 25°C for 6 months

Store at 4°C for 12 months

Store at -20°C for 24 months



Sale Office : No. 16, Sayeh Building, Nafisi Ave., Phase 1, Ekbatan, Tehran, Iran

Tell : +9821 44661950 , +9821 44661903

www.sinaclon.com

www.sinaclon.ir

order@sinaclon.com